
Antibodies Product Datasheet

Product Name:	AKT mouse Monoclonal Antibody(10D6)
Catalog Number:	N1600
Reactivity:	Human,Mouse,Rat,Drosophila
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of AKT at AA range of 400-480
Concentration:	1 mg/ml
Gene Name:	AKT1/AKT2/AKT3
Protein Name:	AKT1
Other Name:	AKT1
Human Swiss Prot No.:	P31749
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AKT protein detects endogenous levels of AKT1
Storage and Stability:	-20°C/1 year