

Antibodies Product Datasheet

Product Name:	NFATc3 Antibody
Catalog Number:	B0522
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NFAT4. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	NFATC3
Protein Name:	Nuclear factor of activated T-cells cytoplasmic 3
Other Name:	NFATC3; NFAT4; Nuclear factor of activated T-cells; cytoplasmic 3; NF-ATc3; NFATc3; NFATx; T-cell transcription factor NFAT4; NF-AT4
Human Swiss Prot No.:	Q12968
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NFATc3 Polyclonal Antibody detects endogenous levels of NFATc3 protein.
Storage and Stability:	-20°C/1 year