

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human Fascin (ABT171) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0340
<b>Clone ID:</b>	ABT171
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC, WB
<b>Application Dilutions</b>	IHC 1:100-200, WB 1:500-1000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	FSCN1
<b>Gene Name:</b>	Q16658
<b>Other Name:</b>	55 kDa actin bundling protein, 55 kDa actin-bundling protein, Actin bundling protein, actin bundling protein, 55-KD, FAN 1, FAN1, Fascin 1, Fascin actin bundling protein 1, Fascin, Fascin homolog 1 actin bundling protein (Strongylocentrotus purpuratus), Fascin homolog 1, Fascin, sea urchin, homolog of, 1, Fascin1, FLJ38511, FSCN 1, FSCN1, FSCN1_HUMAN, HSN, p55, Singed (Drosophila) like (sea urchin fascin homolog like), Singed drosophila homolog like, Singed like (fascin homolog sea urchin) (Drosophila), Singed like (fascin homolog sea urchin), Singed like protein, Singed, drosophila, homolog of, Singed-like protein, SNL, Strongylocentrotus purpuratus
<b>Isotype</b>	IgG1, Kappa
<b>Localization:</b>	Tonsil
<b>Positive Tissue Control:</b>	Cytoplasmic
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human Fascin protein. In western blotting of Hela cell lysate, the antibody can label a 55 kDa band corresponding to Fascin.
<b>Storage and Stability:</b>	-20°C/1 year