

Antibodies Product Datasheet

Product Name:	Tra-2 α Antibody
Catalog Number:	C11000
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TRA-2 alpha. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	TRA2A
Protein Name:	Transformer-2 protein homolog alpha
Other Name:	TRA2A; Transformer-2 protein homolog alpha; TRA-2 alpha; TRA2-alpha; Transformer-2 protein homolog A
Human Swiss Prot No.:	Q13595
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Tra-2 α Polyclonal Antibody detects endogenous levels of Tra-2 α protein.
Storage and Stability:	-20°C/1 year