

Antibodies Product Datasheet

Product Name:	Anti-S100P [ABT-S100P]
Catalog Number:	V0124
Clone ID:	ABT-S100P
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	P25815
Gene Name:	S100P
Other Name:	Protein S100-P; S100 calcium-binding protein P; S100E
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-S100P [ABT-S100P] detects endogenous levels of total S100P protein.
Storage and Stability:	-20°C/1 year