

Antibodies Product Datasheet

Product Name:	CD22 Antibody
Catalog Number:	R12-2069
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	CD22 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD22.
Concentration:	1 mg/ml
Gene Name:	CD22
Other Name:	B-cell receptor CD22; B-lymphocyte cell adhesion molecule; Sialic acid-binding Ig-like lectin 2; T-cell surface antigen Leu-14
Human Swiss Prot No.:	P20273
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CD22 antibody detects endogenous levels of CD22.
Storage and Stability:	Stable at -20°C for at least 1 year.