

Antibodies Product Datasheet

Product Name:	Histone H2B Antibody
Catalog Number:	B0791
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Histone H2B. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	HIST1H2BB
Protein Name:	Histone H2B type 1-B
Other Name:	HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f
Human Swiss Prot No.:	P33778
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Histone H2B Polyclonal Antibody detects endogenous levels of Histone H2B protein.
Storage and Stability:	-20°C/1 year