

Antibodies Product Datasheet

Product Name:	NFκB-p100 Antibody
Catalog Number:	B7164
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IP, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NF-kappaB p100/p52. AA range:833-882
Concentration:	1 mg/ml
Gene Name:	NFKB2
Protein Name:	Nuclear factor NF-kappa-B p100 subunit
Other Name:	NFKB2; LYT10; Nuclear factor NF-kappa-B p100 subunit; DNA-binding factor KBF2; H2TF1; Lymphocyte translocation chromosome 10 protein; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2; Oncogene Lyt-10; Lyt10
Human Swiss Prot No.:	Q00653
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NFκB-p100 Polyclonal Antibody detects endogenous levels of NFκB-p100 protein.
Storage and Stability:	-20°C/1 year