

Antibodies Product Datasheet

Product Name:	CCBP2 Antibody
Catalog Number:	R12-2573
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	CCBP2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCBP2.
Concentration:	1 mg/ml
Gene Name:	CCBP2
Other Name:	Chemokine-binding protein 2; Chemokine-binding protein D6; C-C chemokine receptor D6; Chemokine receptor CCR-9; Chemokine receptor CCR-10; CCBP2; CCR10, CMKBR9;
Human Swiss Prot No.:	O00590
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CCBP2 Antibody antibody detects endogenous levels of CCBP2.
Storage and Stability:	Stable at -20°C for at least 1 year.