

Antibodies Product Datasheet

Product Name:	TRAP100 Antibody
Catalog Number:	C10176
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MED24. AA range:801-850
Concentration:	1 mg/ml
Gene Name:	MED24
Protein Name:	Mediator of RNA polymerase II transcription subunit 24
Other Name:	MED24; ARC100; CRSP4; DRIP100; KIAA0130; THRAP4; TRAP100; Mediator of RNA polymerase II transcription subunit 24; Activator-recruited cofactor 100 kDa component; ARC100; Cofactor required for Sp1 transcriptional activation subunit 4; CRSP c
Human Swiss Prot No.:	O75448
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TRAP100 Polyclonal Antibody detects endogenous levels of TRAP100 protein.
Storage and Stability:	-20°C/1 year