

Antibodies Product Datasheet

Product Name:	ELOVL1 Antibody
Catalog Number:	C15617
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ELOVL1. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	ELOVL1
Protein Name:	Elongation of very long chain fatty acids protein 1
Other Name:	ELOVL1; SSC1; CGI-88; Elongation of very long chain fatty acids protein 1; 3-keto acyl-CoA synthase ELOVL1; ELOVL fatty acid elongase 1; ELOVL FA elongase 1
Human Swiss Prot No.:	Q9BW60
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ELOVL1 Polyclonal Antibody detects endogenous levels of ELOVL1 protein.
Storage and Stability:	-20°C/1 year