

Antibodies Product Datasheet

Product Name:	TMEM132A Antibody
Catalog Number:	C49251
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human TMEM132A
Concentration:	1 mg/ml
Gene Name:	TMEM132A HSPA5BP1 KIAA1583
Protein Name:	Transmembrane protein 132A (HSPA5-binding protein 1)
Other Name:	
Human Swiss Prot No.:	Q24JP5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TMEM132A at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year