

Antibodies Product Datasheet

Product Name:	Anti-Human Gastrin (ABT494) Mouse Monoclonal Antibody
Catalog Number:	V0346
Clone ID:	ABT494
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	GAST
Gene Name:	P01350
Other Name:	Gastric mucin 6, Gastric Mucin, Gastric mucin-6, MUC 6, MUC-6, Muc6, MUC6 Fragment, MUC6 mucin, MUC6 mucin Fragment, MUC6_HUMAN, Mucin 6, Mucin 6 gastric, Mucin 6 oligomeric mucus/gel forming, Mucin glycoprotein Fragment, Mucin-6, Secretory mucin MUC6 Fragment
Isotype	IgG1, Kappa
Localization:	Pyloric region of the stomach
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Gastrin protein.
Storage and Stability:	-20°C/1 year