

---

## Antibodies Product Datasheet

<b>Product Name:</b>	BIOTIN*Polyclonal Goat Anti Human IgG( $\gamma$ chain)
<b>Catalog Number:</b>	SA0301
<b>Reactivity:</b>	HumanIgG( $\gamma$ chain)
<b>Conjugate:</b>	Biotin
<b>Applications:</b>	WB,IHC-p,IF/ICC,ELISA
<b>Dilutions:</b>	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Goat
<b>Purification:</b>	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
<b>Concentration:</b>	4-5mg/mL
<b>Other Name:</b>	BIOTIN*Polyclonal Goat Anti-Human IgG( $\gamma$ chain)
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage and Stability:</b>	-20°C/1 year