

Antibodies Product Datasheet

Product Name:	TSPAN15 Antibody
Catalog Number:	C50321
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human TSPAN15
Concentration:	1 mg/ml
Gene Name:	TSPAN15 NET7 TM4SF15 UNQ677/PRO1311
Protein Name:	Tetraspanin-15 (Tspan-15) (Tetraspan NET-7) (Transmembrane 4 superfamily member 15)
Other Name:	
Human Swiss Prot No.:	O95858
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TSPAN15 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year