

Antibodies Product Datasheet

Product Name:	TNRC4 Antibody
Catalog Number:	C50742
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human TNRC4
Concentration:	1 mg/ml
Gene Name:	CELF3 BRUNOL1 CAGH4 ERDA4 TNRC4
Protein Name:	CUGBP Elav-like family member 3 (CELF-3) (Bruno-like protein 1) (CAG repeat protein 4) (CUG-BP- and ETR-3-like factor 3) (ELAV-type RNA-binding protein 1) (ETR-1) (Expanded repeat domain protein CAG/CTG 4) (RNA-binding protein BRUNOL-1) (Trinucleotide rep
Other Name:	
Human Swiss Prot No.:	Q5SZQ8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TNRC4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year