

Antibodies Product Datasheet

Product Name:	Rpb1 CTD (Phospho-Thr4) Antibody
Catalog Number:	A9751
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Rpb1 CTD (Thr4)
Concentration:	1 mg/ml
Gene Name:	POLR2A POLR2
Protein Name:	Rpb1 CTD (Thr4)
Other Name:	DNA-directed RNA polymerase II subunit RPB1 (RNA polymerase II subunit B1) (EC 2.7.7.6) (DNA-directed RNA polymerase II subunit A) (DNA-directed RNA polymerase III largest subunit) (RNA-directed RNA polymerase II subunit RPB1) (EC 2.7.7.48)
Human Swiss Prot No.:	P24928
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat Rpb1 CTD (phospho-Thr4)
Storage and Stability:	-20°C/1 year