

Antibodies Product Datasheet

Product Name:	NFκB-p65 (Acetyl-Lys310) Antibody
Catalog Number:	D0018
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	WB 1:500-2000, IHC-p 1:50-300, IF 1:50-300
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NF-κappaB p65 around the acetylated site of Lys310. AA range:275-324
Concentration:	1 mg/ml
Gene Name:	RELA
Protein Name:	Transcription factor p65
Other Name:	RELA; NFKB3; Transcription factor p65; Nuclear factor NF-κappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3
Human Swiss Prot No.:	Q04206
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-NFκB-p65 (K310) Polyclonal Antibody detects endogenous levels of NFκB-p65 protein only when acetylated at K310.
Storage and Stability:	-20°C/1 year