

Antibodies Product Datasheet

Product Name:	FA5 (heavy chain, Cleaved-Arg737) Antibody
Catalog Number:	L0426
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FA5 (heavy chain, Cleaved-Arg737)
Concentration:	1 mg/ml
Gene Name:	F5
Protein Name:	FA5 (heavy chain, Cleaved-Arg737)
Other Name:	Coagulation factor V (Activated protein C cofactor) (Proaccelerin, labile factor) [Cleaved into: Coagulation factor V heavy chain; Coagulation factor V light chain]
Human Swiss Prot No.:	P12259
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat FA5 (heavy chain, Cleaved-Arg737, protein was cleaved amino acid sequence between 737-738)
Storage and Stability:	-20°C/1 year