

Antibodies Product Datasheet

Product Name:	GPX5 Antibody
Catalog Number:	C48344
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GPX5
Concentration:	1 mg/ml
Gene Name:	GPX5
Protein Name:	Epididymal secretory glutathione peroxidase (EC 1.11.1.9) (Epididymis-specific glutathione peroxidase-like protein) (EGLP) (Glutathione peroxidase 5) (GPx-5) (GSHPx-5)
Other Name:	
Human Swiss Prot No.:	O75715
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GPX5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year