
Antibodies Product Datasheet

Product Name:	COT Antibody
Catalog Number:	B0064
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human COT. AA range:256-305
Concentration:	1 mg/ml
Gene Name:	MAP3K8
Protein Name:	Mitogen-activated protein kinase kinase kinase 8
Other Name:	MAP3K8; COT; ESTF; Mitogen-activated protein kinase kinase kinase 8; Cancer Osaka thyroid oncogene; Proto-oncogene c-Cot; Serine/threonine-protein kinase cot; Tumor progression locus 2; TPL-2
Human Swiss Prot No.:	P41279
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	COT Polyclonal Antibody detects endogenous levels of COT protein.
Storage and Stability:	-20°C/1 year