

## Antibodies Product Datasheet

<b>Product Name:</b>	RBL2 (Phospho Thr642) Antibody
<b>Catalog Number:</b>	C51127
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC, WB
<b>Dilutions:</b>	WB 1:500-2000 IHC 1:50-200
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Rabbit, polyclonal
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human RBL2 (Phospho Thr642)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	RBL2 RB2
<b>Protein Name:</b>	Retinoblastoma-like protein 2 (130 kDa retinoblastoma-associated protein) (p130) (Retinoblastoma-related protein 2) (RBR-2) (pRb2)
<b>Other Name:</b>	Retinoblastoma-like protein 2 (130 kDa retinoblastoma-associated protein) (p130) (Retinoblastoma-related protein 2) (RBR-2) (pRb2)
<b>Human Swiss Prot No.:</b>	Q08999
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of RBL2 (Phospho Thr642) Rabbit pAb at Human, Mouse, Rat
<b>Storage and Stability:</b>	Store at -20° C/1 year