

---

## Antibodies Product Datasheet

|                               |  |
|-------------------------------|--|
| <b>Product Name:</b>          | SGK1 mouse Monoclonal Antibody(4D12)   |
| <b>Catalog Number:</b>        | N1597  |
| <b>Reactivity:</b>            | Human,Rat,Mouse  |
| <b>Applications:</b>          | IHC-p,IF(paraffin section)   |
| <b>Dilutions:</b>             | IHC 1:100-200  |
| <b>Modification:</b>          | Unmodified/ Total  |
| <b>Source:</b>                | Monoclonal Mouse   |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| <b>Immunogen:</b>             | Synthetic Peptide of SGK1 at AA range of 350-430   |
| <b>Concentration:</b>         | 1 mg/ml  |
| <b>Gene Name:</b>             | SGK1   |
| <b>Protein Name:</b>          | SGK1   |
| <b>Other Name:</b>            | SGK1   |
| <b>Human Swiss Prot No.:</b>  | O00141   |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                    |
| <b>Specificity:</b>           | SGK1 protein detects endogenous levels of SGK1   |
| <b>Storage and Stability:</b> | -20°C/1 year   |