

Antibodies Product Datasheet

Product Name:	FH Antibody
Catalog Number:	C47447
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Fumarate hydratase, mitochondrial
Concentration:	
Gene Name:	FH
Protein Name:	Fumarate hydratase, mitochondrial
Other Name:	Fumarate hydratase, mitochondrial (Fumarase) (EC 4.2.1.2); YM3222
Human Swiss Prot No.:	P07954
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous fumarase proteins.
Storage and Stability:	-20°C/1 year