
Antibodies Product Datasheet

Product Name:	Insulin mouse mAb(PT0156)
Catalog Number:	N1740
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500, IF 1:100-500
Modification:	Unmodified/ Total
Source:	Mouse,Monoclonal
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Insulin
Concentration:	
Gene Name:	INS
Protein Name:	Insulin
Other Name:	Insulin [Cleaved into: Insulin B chain; Insulin A chain]
Human Swiss Prot No.:	P01308
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.26% sodium azide.
Specificity:	This antibody detects endogenous levels of Insulin at Human
Storage and Stability:	-20°C/1 year