

Antibodies Product Datasheet

Product Name:	GNPDA2 Antibody
Catalog Number:	C48932
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GNPDA2
Concentration:	1 mg/ml
Gene Name:	GNPDA2 GNP2
Protein Name:	Glucosamine-6-phosphate isomerase 2 (EC 3.5.99.6) (Glucosamine-6-phosphate deaminase 2) (GNPDA 2) (GlcN6P deaminase 2) (Glucosamine-6-phosphate isomerase SB52)
Other Name:	
Human Swiss Prot No.:	Q8TDQ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GNPDA2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year