
Antibodies Product Datasheet

Product Name:	TAT Antibody
Catalog Number:	R12-2376
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAT. AA range:255-304
Concentration:	1 mg/ml
Gene Name:	TAT
Protein Name:	Tyrosine aminotransferase
Other Name:	TAT; Tyrosine aminotransferase; TAT; L-tyrosine:2-oxoglutarate aminotransferase
Human Swiss Prot No.:	P17735
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TAT Polyclonal Antibody detects endogenous levels of TAT protein.
Storage and Stability:	-20°C/1 year