

Antibodies Product Datasheet

Product Name:	Anti-Human Stathmin (ABT-Stathmin) Mouse Monoclonal Antibody
Catalog Number:	V0159
Clone ID:	ABT-Stathmin
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	STMN1
Gene Name:	P16949
Other Name:	C1orf215, Lag, LAP 18, LAP18, Leukemia associated phosphoprotein p18, Leukemia-associated phosphoprotein p18, Metablastin, Oncoprotein 18, OP 18, Op18, p18, p19, Phosphoprotein 19, Phosphoprotein p19, pp17, pp19, PR22, Pr22 protein, Prosolin, Protein Pr22, SMN, Stathmin, Stathmin1, STMN 1, Stmn1, STMN1_HUMAN
Isotype	IgG2b, Kappa
Localization:	Tonsil
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Stathmin protein. In western blotting of Hela, HepG2 and HEK293 cell lysates, the antibody can label a 17 kDa band corresponding to Stathmin.
Storage and Stability:	-20°C/1 year