

Antibodies Product Datasheet

Product Name:	CK8 Monoclonal Antibody(8G8)
Catalog Number:	N1320
Reactivity:	Human,Mouse,Rat
Applications:	IF/ICC,WB,IHC-p
Dilutions:	IF: 1:50-200 WB: 1:2000-5000 IHC: 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CK8
Concentration:	
Gene Name:	KRT8
Protein Name:	Keratin type II cytoskeletal 8
Other Name:	KRT8; CYK8; Keratin, type II cytoskeletal 8; Cytokeratin-8; CK-8; Keratin-8; K8; Type-II keratin Kb8
Human Swiss Prot No.:	P05787
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CK8 proteins.
Storage and Stability:	-20°C/1 year