

Antibodies Product Datasheet

Product Name:	VE-Cadherin mouse Monoclonal Antibody(3G8)
Catalog Number:	N1705
Reactivity:	Human,Mouse,Rat
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of VE-Cadherin
Concentration:	1 mg/ml
Gene Name:	CDH5
Protein Name:	Cadherin-5 (7B4 antigen) (Vascular endothelial cadherin) (VE-cadherin) (CD antigen CD144)
Other Name:	Cadherin-5 (7B4 antigen) (Vascular endothelial cadherin) (VE-cadherin) (CD antigen CD144)
Human Swiss Prot No.:	P33151
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous VE-Cadherin protein
Storage and Stability:	-20°C/1 year