
Antibodies Product Datasheet

Product Name:	Caldesmon Antibody
Catalog Number:	R12-2566
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Caldesmon Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Caldesmon.
Concentration:	1 mg/ml
Gene Name:	CALD1
Other Name:	CAD; CALD1; CDM; L-caldesmon; Non-muscle caldesmon
Human Swiss Prot No.:	Q05682
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Caldesmon Antibody antibody detects endogenous levels of Caldesmon.
Storage and Stability:	Stable at -20°C for at least 1 year.