
Antibodies Product Datasheet

Product Name:	CDC25C Monoclonal Antibody
Catalog Number:	N0141
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CDC25C expressed in E. Coli.
Concentration:	
Gene Name:	CDC25C
Protein Name:	M-phase inducer phosphatase 3
Other Name:	CDC25C; M-phase inducer phosphatase 3; Dual specificity phosphatase Cdc25C
Human Swiss Prot No.:	P30307
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CDC25C Monoclonal Antibody detects endogenous levels of CDC25C protein.
Storage and Stability:	-20°C/1 year