

Antibodies Product Datasheet

Product Name:	RYBP Antibody
Catalog Number:	C50541
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RYBP
Concentration:	1 mg/ml
Gene Name:	RYBP DEDAF YEAF1
Protein Name:	RING1 and YY1-binding protein (Apoptin-associating protein 1) (APAP-1) (Death effector domain-associated factor) (DED-associated factor) (YY1 and E4TF1-associated factor 1)
Other Name:	
Human Swiss Prot No.:	Q8N488
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RYBP at Human, Mouse
Storage and Stability:	Store at -20° C/1 year