

## Antibodies Product Datasheet

<b>Product Name:</b>	Olfactory receptor 10G9 Antibody
<b>Catalog Number:</b>	G420
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	IF/ICC, ELISA
<b>Dilutions:</b>	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human OR10G9. AA range:231-280
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	OR10G9
<b>Protein Name:</b>	Olfactory receptor 10G9
<b>Other Name:</b>	OR10G9; OR10G10P; Olfactory receptor 10G9; Olfactory receptor 10G10
<b>Human Swiss Prot No.:</b>	Q8NGN4
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Olfactory receptor 10G9 Polyclonal Antibody detects endogenous levels of Olfactory receptor 10G9 protein.
<b>Storage and Stability:</b>	-20°C/1 year