

Antibodies Product Datasheet

Product Name:	DCAMKL2 Antibody
Catalog Number:	C11345
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DCLK2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	DCLK2
Protein Name:	Serine/threonine-protein kinase DCLK2
Other Name:	DCLK2; DCAMKL2; DCDC3B; DCK2; Serine/threonine-protein kinase DCLK2; CaMK-like CREB regulatory kinase 2; CL2; CLICK-II; CLICK2; Doublecortin domain-containing protein 3B; Doublecortin-like and CAM kinase-like 2; Doublecortin-like kinase 2
Human Swiss Prot No.:	Q8N568
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DCAMKL2 Polyclonal Antibody detects endogenous levels of DCAMKL2 protein.
Storage and Stability:	-20°C/1 year