

Antibodies Product Datasheet

Product Name:	B3GNT3 Antibody
Catalog Number:	C48462
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human B3GNT3
Concentration:	1 mg/ml
Gene Name:	B3GNT3 B3GALT8 TMEM3 UNQ637/PRO1266
Protein Name:	UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase 3 (BGnT-3) (Beta-1,3-Gn-T3) (Beta-1,3-N-acetylglucosaminyltransferase 3) (Beta3Gn-T3) (EC 2.4.1.-) (Beta-1,3-galactosyltransferase 8) (Beta-1,3-GalTase 8) (Beta3Gal-T8) (Beta3GalT8) (b3Gal-T8) (B
Other Name:	
Human Swiss Prot No.:	Q9Y2A9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of B3GNT3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year