

Antibodies Product Datasheet

Product Name:	SRF (Phospho-Thr159) Antibody
Catalog Number:	A8427
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SRF around the phosphorylation site of Thr159. AA range:126-175
Concentration:	1 mg/ml
Gene Name:	SRF
Protein Name:	Serum response factor
Other Name:	SRF; Serum response factor; SRF
Human Swiss Prot No.:	P11831
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-SRF (T159) Polyclonal Antibody detects endogenous levels of SRF protein only when phosphorylated at T159.
Storage and Stability:	-20°C/1 year