

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human Perforin (ABT355) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0434
<b>Clone ID:</b>	ABT355
<b>Reactivity:</b>	Human
<b>Applications:</b>	IHC, WB
<b>Application Dilutions</b>	IHC 1:100-200, WB 1:500-1000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	PRF1
<b>Gene Name:</b>	P14222
<b>Other Name:</b>	Cytolysin, FLH2, HPLH2, Lymphocyte pore-forming protein, P1, PERF_HUMAN, perforin 1 (pore forming protein), Perforin 1, Perforin-1, PFP, PGFL, PIGF, PIGF-2, PLGF, Pore forming protein, prf1, SHGC-10760
<b>Isotype</b>	IgG1, Kappa
<b>Localization:</b>	Spleen
<b>Positive Tissue Control:</b>	Cytoplasmic, Membranous
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human Perforin protein.
<b>Storage and Stability:</b>	-20°C/1 year