

Antibodies Product Datasheet

Product Name:	PIP4K2A Antibody
Catalog Number:	C49616
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PIP4K2A
Concentration:	1 mg/ml
Gene Name:	PIP4K2A PIP5K2 PIP5K2A
Protein Name:	Phosphatidylinositol 5-phosphate 4-kinase type-2 alpha (EC 2.7.1.149) (1-phosphatidylinositol 5-phosphate 4-kinase 2-alpha) (Diphosphoinositide kinase 2-alpha) (PIP5KIII) (Phosphatidylinositol 5-phosphate 4-kinase type II alpha) (PI(5)P 4-kinase type II a
Other Name:	
Human Swiss Prot No.:	P48426
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PIP4K2A at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year