

Antibodies Product Datasheet

Product Name:	MARCH5 Antibody
Catalog Number:	C15550
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MARCH5. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	MARCH5
Protein Name:	E3 ubiquitin-protein ligase MARCH5
Other Name:	MARCH5; RNF153; E3 ubiquitin-protein ligase MARCH5; Membrane-associated RING finger protein 5; Membrane-associated RING-CH protein V; MARCH-V; Mitochondrial ubiquitin ligase; MITOL; RING finger protein 153
Human Swiss Prot No.:	Q9NX47
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MARCH5 Polyclonal Antibody detects endogenous levels of 40607 protein.
Storage and Stability:	-20°C/1 year