

Antibodies Product Datasheet

Product Name:	MyoD (Acetyl-Lys99/K102) Antibody
Catalog Number:	D12218
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MyoD (Acetyl Lys99/K102)
Concentration:	1 mg/ml
Gene Name:	MYOD1 BHLHC1 MYF3 MYOD
Protein Name:	MyoD (Acetyl Lys99/K102)
Other Name:	Myoblast determination protein 1 (Class C basic helix-loop-helix protein 1) (bHLHc1) (Myogenic factor 3) (Myf-3)
Human Swiss Prot No.:	P15172
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human MyoD (Acetyl Lys99/K102)
Storage and Stability:	-20°C/1 year