

## Antibodies Product Datasheet

<b>Product Name:</b>	Smad3 Monoclonal Antibody(6F7)
<b>Catalog Number:</b>	N1560
<b>Reactivity:</b>	Human,Rat,Mouse
<b>Applications:</b>	IHC-p,IF(paraffin section)
<b>Dilutions:</b>	IHC 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Recombinant Protein of Smad3
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	SMAD3
<b>Protein Name:</b>	Mothers against decapentaplegic homolog 3 (MAD homolog 3) (Mad3) (Mothers against DPP homolog 3) (hMAD-3) (JV15-2) (SMAD family member 3) (SMAD 3) (Smad3) (hSMAD3)
<b>Other Name:</b>	SMAD3; MADH3; Mothers against decapentaplegic homolog 3; MAD homolog 3; Mad3; Mothers against DPP homolog 3; hMAD-3; JV15-2; SMAD family member 3; SMAD 3; Smad3; hSMAD3
<b>Human Swiss Prot No.:</b>	P84022
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Smad3 protein detects endogenous levels of Smad3
<b>Storage and Stability:</b>	-20°C/1 year