

Antibodies Product Datasheet

Product Name:	Fractalkine Receptor Antibody
Catalog Number:	C20458
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Fractalkine Receptor . at AA range: 120-200
Concentration:	1 mg/ml
Gene Name:	CX3CR1
Protein Name:	CX3C chemokine receptor 1
Other Name:	CX3CR1; CMKBRL1; GPR13; CX3C chemokine receptor 1; C-X3-C CKR-1; CX3CR1; Beta chemokine receptor-like 1; CMK-BRL-1; CMK- BRL1; Fractalkine receptor; G-protein coupled receptor 13; V28
Human Swiss Prot No.:	P49238
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fractalkine Receptor Polyclonal Antibody detects endogenous levels of Fractalkine Receptor protein.
Storage and Stability:	-20°C/1 year