

---

## Antibodies Product Datasheet

<b>Product Name:</b>	BCL-2 Mouse Monoclonal Antibody(2F1)
<b>Catalog Number:</b>	N1721
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB,IHC
<b>Dilutions:</b>	IHC-p 1:50-200, WB 1:500-2000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human BCL-2
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	BCL2
<b>Protein Name:</b>	BCL-2
<b>Other Name:</b>	Apoptosis regulator Bcl-2
<b>Human Swiss Prot No.:</b>	P10415
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of human BCL-2
<b>Storage and Stability:</b>	-20°C/1 year