

Antibodies Product Datasheet

Product Name:	DNASE1L3 Antibody
Catalog Number:	C50190
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DNASE1L3
Concentration:	1 mg/ml
Gene Name:	DNASE1L3 DHP2 DNAS1L3
Protein Name:	Deoxyribonuclease gamma (DNase gamma) (EC 3.1.21.-) (DNase I homolog protein DHP2) (Deoxyribonuclease I-like 3) (DNase I-like 3) (Liver and spleen DNase) (LS-DNase) (LSD)
Other Name:	
Human Swiss Prot No.:	Q13609
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DNASE1L3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year