

Antibodies Product Datasheet

Product Name:	NPY5-R Antibody
Catalog Number:	G412
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NPY5R. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	NPY5R
Protein Name:	Neuropeptide Y receptor type 5
Other Name:	NPY5R; NPYR5; Neuropeptide Y receptor type 5; NPY5-R; NPY-Y5 receptor; NPYY5-R; Y5 receptor
Human Swiss Prot No.:	Q15761
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NPY5-R Polyclonal Antibody detects endogenous levels of NPY5-R protein.
Storage and Stability:	-20°C/1 year