
Antibodies Product Datasheet

Product Name:	TIMP-4 Antibody
Catalog Number:	C30897
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from TIMP-4 at AA range: 174-223
Concentration:	1 mg/ml
Gene Name:	TIMP4
Protein Name:	TIMP-4
Other Name:	Metalloproteinase inhibitor 4 (Tissue inhibitor of metalloproteinases 4) (TIMP-4)
Human Swiss Prot No.:	Q99727
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TIMP-4 Polyclonal Antibody detects endogenous levels of TIMP-4
Storage and Stability:	-20°C/1 year