
Antibodies Product Datasheet

Product Name:	Cleaved-MMP-27 (Y99) Antibody
Catalog Number:	L0334
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MMP27. AA range:80-129
Concentration:	1 mg/ml
Gene Name:	MMP27
Protein Name:	Matrix metalloproteinase-27
Other Name:	MMP27; Matrix metalloproteinase-27; MMP-27
Human Swiss Prot No.:	Q9H306
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-MMP-27 (Y99) Polyclonal Antibody detects endogenous levels of fragment of activated MMP-27 protein resulting from cleavage adjacent to Y99.
Storage and Stability:	-20°C/1 year