
Antibodies Product Datasheet

Product Name:	JAK2 mouse Monoclonal Antibody(9B4)
Catalog Number:	N1689
Reactivity:	Human,Mouse,Rat
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of JAK2
Concentration:	1 mg/ml
Gene Name:	JAK2
Protein Name:	Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2)
Other Name:	Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2)
Human Swiss Prot No.:	O60674
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous JAK2 protein
Storage and Stability:	-20°C/1 year