

Antibodies Product Datasheet

Product Name:	SLC37A4 Antibody
Catalog Number:	C50736
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC37A4
Concentration:	1 mg/ml
Gene Name:	SLC37A4 G6PT G6PT1 PRO0685 TRG19
Protein Name:	Glucose-6-phosphate translocase (Glucose-5-phosphate transporter) (Solute carrier family 37 member 4) (Transformation-related gene 19 protein) (TRG-19)
Other Name:	
Human Swiss Prot No.:	O43826
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC37A4 at Human
Storage and Stability:	Store at -20° C/1 year